

Introducing the **ProKase™** for
Highly Portable Flood Source
Storage and Handling.



ProKase™
FLOOD SHIELD



Eckert & Ziegler
Isotope Products



More... Mobile Stable Access




Wide Clamshell Opening
for Easy Access



Ergonomic Top Handle
for Easier Rolling



Pull Handle
to Exit Confined Spaces



Unique Flared Base
for Improved Stability



Wheel Brake
for Added Safety



Four Swivel Wheels
for Greater Mobility





ProKem Technology™ = Better Quality Control

ProKem Technology begins with our proprietary chemistry which results in a unique purification and production process that produces cleaner Cobalt for a more uniform source. This yields superior camera calibration, less variation and, ultimately, better quality control.



Eckert & Ziegler
Isotope Products

Eckert & Ziegler Isotope Products
24937 Avenue Tibbitts
Valencia, CA 91355

USA
Phone: +1 (661) 309-1010
Fax: +1 (661) 257-8303
Toll Free: +1 (866) 476-9767

nucmedsales@ezag.com
www.ezag.com

Eckert & Ziegler Isotope Products GmbH
Betriebstätte Braunschweig
Gieselweg 1
38110 Braunschweig

Germany
Phone: +49 5307 932-0
Fax: +49 5307 932-272

ipe@ezag.de
www.ezag.com

Eckert & Ziegler Isotope Products Sarl
12 Avenue des Tropiques
Hightec Sud - Bâtiment B
91955 Courtaboeuf Cedex

France
Phone: +33 1 64 86 22 22
Fax: +33 1 69 86 10 65

jeanrene.caramelle@ezag-france.fr
www.ezag.com